



Falcon R&D Program for
Dermato-Oncology



Identification of patients at high risk for relapse using the Merlin Assay (CP-GEP) in an independent cohort of 930 patients (pts) with stage I/II melanoma who did not undergo sentinel lymph node biopsy.

T. Amaral¹, E. Chatziioannou¹, A. Nuebling¹, T. Sinnberg¹, H. Niessner¹, U. Leiter¹, J. Dwarkasing², T. Arentsen², C. Groß², A. Eggermont^{3,4}, L. Flatz¹, S. Forchhammer¹

¹Skin Cancer Clinical Trials Center, Tuebingen, Germany, ²SkylineDX, Rotterdam, Netherlands, ³University Medical Center Utrecht & Princess Maxima Center, Utrecht, Netherlands, ⁴Comprehensive Cancer Center Munich of the Technical University Munich and the Ludwig Maximilians University, Munich, Germany

Introduction

Sentinel lymph node biopsy (SLNB) is the gold standard for nodal assessment in staging cutaneous melanoma (CM) according to AJCC v8 guideline. More than 80% of pts are negative for nodal metastasis, but most pts who relapse or die from melanoma are initially diagnosed as 'low risk' early-stage. Previously we showed that the clinicopathological-gene expression profiling (CP-GEP) model can stratify stage I-II pts and pts who did not undergo SLNB in low and high-risk of recurrence (Amaral et al., EJC, 2023). Here we investigate CP-GEP ability to stratify pts who did not undergo SLNB for their risk of recurrence in substantial cohort.

Methods

We analysed formalin-fixed paraffin-embedded primary tumor samples of 930 pts with stage I/II CM diagnosed between 2000-2017, included in the Central Malignant Melanoma Registry, who did not

This document is intended solely for informational and educational purposes as a scientific resource and is not intended for commercial use.

receive SLNB. Tumors were analysed blinded to clinical outcome. The CP-GEP model used combines the expression of 8 genes (SERPINE2, GDF15, ITGB3, CXCL8, LOXL4, TGFBR1, PLAT and MLANA) by quantitative reverse transcription polymerase chain reaction with age and Breslow thickness to obtain a binary output: CP-GEP Low- or High-Risk. Relapse-free survival (RFS), distant metastasis free survival (DMFS) and Melanoma Specific Survival (MSS) were evaluated using Kaplan-Meier curves.

Results

We included 930 pts (stage IA-IIc). 41% were females, median age was 64-year-old, median Breslow thickness was 0.5 mm, the majority were not ulcerated (94%). For all pts, the 5-year RFS was 90.9%; 5-year DMFS was 96.9 and 5-year MSS was 97.5%. Median follow-up time was 55 months (RFS). CP-GEP identified 879 pts as Low-Risk and 51 pts as High-Risk. The 5-year RFS rate was 94.6% for CP-GEP Low-Risk pts versus 26.6% for CP-GEP High-Risk patients (HR 25.08; $p < 0.001$). 5-year DMFS was 98.6% vs 62.1% (HR 35.39; $p < 0.001$) for CP-GEP Low-Risk and High-Risk pts, respectively. The 5-year MSS was 99.4% for Low-Risk and 61.7% for High-Risk pts (HR 71.05; $p < 0.001$), capturing 12 out of 16 melanoma specific deaths in the CP-GEP High-Risk group.

Conclusion

This comprehensive study shows that CP-GEP has the potential to stratify pts with early-stage melanoma who did not undergo SLNB based on their risk of recurrence. Pts with CP-GEP Low-Risk have a good long-term survival while pts with CP-GEP High-Risk have a high risk of recurrence. CP-GEP may have the potential to stratify pts beyond SLNB.

Abstract at EADO 2025

This document is intended solely for informational and educational purposes as a scientific resource and is not intended for commercial use.

